



### RUGGEDIZED TABLET OR LAPTOP

- Large 17" Medically Certified Sunlight Readable Tablet with Capacitive Touchscreen
- Lightweight at 4.4 lbs (Tablet) and 4.34 lbs (Laptop) with long battery life (6 hours of continuous use)
- 4G, 5G Enabled for real-time Remote Support & Study Upload to ConnectVET Cloud anytime / anywhere
- Protected by an IP65 Ingress Rating for Dust & Water and Externally with a Rubberized Protective Bumper

### LIGHTWEIGHT RUGGED DETECTOR

- Unbreakable Flexible Substrate Glass-Free TFT Panel
- 99 micron Pixel Pitch (> 5.0 lp/mm resolution)
- IP67 Water & Dust Protection
- Robust Design: Drop Tested at 1.5m
- Up to 8 hours autonomy
- Lightweight only 3.74lbs
- Polycarbonate Protection with Interchangeable handle





### REMOTE IMAGE PREVIEW

- Shows **actual image** being taken, no need to go back to your laptop
- Convenient retake button
- Mounted on our ConnectRAY portable generator

### CONTACT US

 [www.imv-technologies.us](http://www.imv-technologies.us)  
 [info@imvusa.com](mailto:info@imvusa.com)  
 (507) 529-8200  
 (800) 210-9665

### FOLLOW US

 [www.facebook.com/imvimaging.northamerica](https://www.facebook.com/imvimaging.northamerica)  
 [www.instagram.com/imvimaging.na](https://www.instagram.com/imvimaging.na)



ConnectVET is your daily multi-modality tool to share & visualize your medical images but also has the capacity for videos and photos. It was built to communicate with other Veterinarians, clinics, clients, and PACS systems. Truly the single source & affordable Cloud PACS solution for any & every practice!

### MOBILE




Browse, diagnose and mail your studies right from the **palm of your hand**. Available on iOS and Android

### REPORTS



Split-view enhances the creation of **personalized** reports with selected images, checkboxes, and tables while reviewing studies

### LOCAL + CLOUD



Supports **multiple modalities** (CT, US, MRI, Endo...) both offline and online

### AUTO-MODE



Compare L/R-images automatically for an **easier** reviewing process



### CONTACT US

🏠 [www.imv-technologies.us](http://www.imv-technologies.us)  
 ✉ [info@imvusa.com](mailto:info@imvusa.com)  
 ☎ (507) 529-8200  
 (800) 210-9665

### FOLLOW US

📘 [www.facebook.com/imvimaging.northamerica](https://www.facebook.com/imvimaging.northamerica)  
 📷 [www.instagram.com/imvimaging.na](https://www.instagram.com/imvimaging.na)

